

2023

NURSE CALL SYSTEMS CATALOGUE



www.bellme.co.za

Providing Wireless Nurse Call Systems to South Africa Since 2012.
Allowing instant attention of the nurses with the simple press of a button.

TABLE OF CONTENTS

<u>INTRODUCTION</u>	03.
<u>OUR VALUES</u>	04.
<u>PRODUCT OVERVIEW</u>	05.
<u>HOW BELLME WORKS</u>	06.
<u>APPLICATIONS</u>	07.
<u>SINGLE KEY CALL BUTTON</u>	08.
<u>SINGLE KEY CANCEL BUTTON</u>	11.
<u>CALL & CANCEL WALL UNIT</u>	12.
<u>CALL & CANCEL WALL UNIT WITH CORD</u>	14.
<u>LED LIGHT RECEIVER</u>	16.
<u>WRISTWATCH RECEIVER</u>	18.
<u>B-99P DISPLAY PANEL</u>	20.
<u>B-99E DISPLAY PANEL</u>	22.
<u>B-99G DISPLAY PANEL</u>	24.
<u>B-901 DISPLAY PANEL</u>	26.
<u>W900 DESKTOP SOFTWARE</u>	28.
<u>SIGNAL REPEATER</u>	30.
<u>USB RECEIVER & SOFTWARE</u>	32.
<u>TESTIMONIALS</u>	34.
<u>CLIENTS</u>	37.
<u>CONTACT US</u>	38.

INTRODUCTION

Providing Wireless Nurse Call Systems that are fast becoming the answer to many hospitals and frail care centres in SA, BellMe was established in 2012 by Ushir Shah in his search for a solution to simplify the way various industries receive notifications of when their customer or clients require assistance.

His aim was to increase the reaction time of staff in the case of emergencies by introducing a wireless nurse call system to retirement homes and hospitals.

ABOUT US



Villa Sunfield in Durban North was BellMe's very first nurse call systems client and to this day, they are still making use of the same system. Since then, we have supplied many big clients such as Clicks, Dischem, Akeso Psychiatric Hospitals and many more.

OUR VALUES



Happiness:

Ensuring our team is always happy, motivated and inspired to work, creates an amazing energy and thus ensuring that energy is passed onto you.

Quality:

Ensuring you receive the best value in terms of product and service by making sure all agents and staff are properly trained and all equipment is thoroughly tested before dispatched.

Accountability:

Our extensive warranty indicates trust in our product. We call on you to hold us accountable for this just as we do ourselves and our team for providing excellent customer service.

PRODUCT OVERVIEW

Imagine the possibility of improved response time, minimal cables and instant notification of an emergency with the push of a button. With BellMe we can offer you all of this and so much more.

Quick Installation, Easy to Use:
Self-installation is a breeze. Simply mount the buttons on a flat surface, plug in the display and you are good to go.


Agents Available in Your Area:
Enquire with us for more information regarding agents in your area to conduct an on-site inspection and quotation.



No Hidden Costs:
Your Nurse Call System is a once-off purchase and comes with the basic equipment required for it to run smoothly. No monthly fees.



HOW BELLME WORKS




CALL BUTTON

The notification is wirelessly transmitted from the CALL BUTTON to:



DISPLAY PANEL

The nurse or matron via a DISPLAY PANEL at the nurse's station



WRISTWATCH

A WRISTWATCH gives off an alert sound and vibration in case the nurse is away from the station during a medical emergency.



LED LIGHTS

On the way to the patient or resident, the nurse may also quickly recognize which room to go into via LED LIGHTS flashing above the door.



CANCEL BUTTON

Clear the call from the display panel and the LED Light by pressing the CANCEL BUTTON.





HOSPITALS

Each patient will have a Call Button mounted by their bedside and in the bathrooms to alert the nurses at the nurse's station that assistance is needed and from where.



FRAIL CARE

Each resident or patient will have a Call Button mounted by their bedside or worn on a lanyard to notify the nurse or caregiver when assistance is required and from where.



RETIREMENT

Residents will have a Call Button that they can wear on lanyard and carry with them to notify the nurse or caregiver when assistance is needed and from who.



CLINICS

Each room will have a Call Button for patients or doctors to press which will notify on a display panel indicating which room the call comes from and where assistance is needed.

SINGLE KEY CALL BUTTON

This Call Button Saves Lives.

A simple push of the button will send out alerts to multiple receivers such as:

- A Display Panel,
- A Wristwatch receiver and
- An LED Receiver.



The button is durable, waterproof and long lasting, ensuring your patient or resident has peace of mind for a lifetime.

BENEFITS

5-YEAR WARRANTY:

We have complete faith that our Buttons will continue to Save Lives.

This is why we recently increased our Warranty on our single key Call Buttons from a 2-year to a 5-year warranty as of May 2018. If your Button happens to give you any issues that are covered under your warranty we will even give you a replacement unit until we resolve your technical issue.

100% WATERPROOF:

This call Button can be taken into high-risk areas such as bathrooms, toilets and showers.

It allows the patient or resident to take their trusty Call Button into the shower or out into the garden so they don't have to fear falling down and waiting for hours before someone finds them.

EASY TO PRESS:

- Large surface area that can be triggered on a soft rubber with a slight press.
- The rubber surface is indented to prevent an accidental press.



DUAL PURPOSE APPLICATION:

1 - Wearable:

A lanyard can be connected to the back-mounting section and can be used around the neck by the patient or resident so that they have the freedom to move around and not be confined to the area that their button is located. They can press the call button wherever they are.

2 - Mountable:

Alternatively, it can be wall mounted as it comes with a wall-mounting bracket. You can drill it into the wall with the mounting bracket or use very strong 2-way tape. The Call Button can be mounted on any flat surface within the room and other high-risk areas such as toilets, bathrooms and showers as they are completely waterproof.



A. Insert the key to release the backing.



B. The backing is mounted to a surface.



C. Slide securely onto the mounted bracket.

FEATURES

LONG LASTING BATTERY USAGE:

- Requires an A23 battery.
- Battery can last up to 1-2 years in a Nursing environment.
- Supplied with a special key to remove the backing and unscrew the casing to change the battery.

PRE-PROGRAMMED WITH A UNIQUE CODE:

- Each call button has a unique code
- The code can be saved to display what you would like.

EASY BATTERY CHANGE

Easily change the A23 battery within three simple steps.



A. Insert the key to release the backing.



B. The backing is mounted to a surface.



C. Remove the 3 screws to replace the battery.



Wear the button on a lanyard around the neck.



Mount the button on any flat surface by the bedside.



Mount the button on any flat surface in bathrooms.

B-A1-WL & B-A1-BC Technical Information	Dimensions: 60x60mm	Radio Frequency: 433MHz
	Battery: A23	Range: 100 - 150m depending on the surroundings

CANCEL BUTTON

MONITOR RESPONSE TIME:

If you would like to have the response time of the Nurses captured via the BellMe Software. The nurse would have to actively cancel the call button either on the screen or a cancel button at the bedside or call point for the call to be removed from the display.



MOUNTABLE:

The Cancel button is mountable exactly like the Call Button. It comes with a mounting bracket and can be mounted with strong 2-way tape or drilled.

CANCEL ONE CALL OR CANCEL ALL:

Sharing a Cancel Button:

This button is usually placed at the door to cancel calls coming from that section/ward. If you are sharing a cancel button with multiple call buttons then your display screen will have to display the same calling information for the multiple call buttons.

Individual Cancel Buttons:

The button can be wall mounted near the Bedside Call Buttons or in the Toilets or Shower.

LONG LASTING BATTERY USAGE:

- Requires an A23 battery.
- Battery can last up to 1-2 years in a Nursing environment.
- Supplied with a special key to remove the backing and unscrew the casing to change the battery.



Call and Cancel Button mounted on the patient or resident's bedside.



Call Button mounted in a Toilet. The Cancel Button is usually mounted at the Toilet Entrance.

CALL AND CANCEL WALL UNIT

A simple push of the button will send out alerts to multiple receivers such as:

- A Display Panel,
- A Wristwatch receiver and
- An LED Receiver.



The large Call Button can send out alerts to multiple receivers such as a Display Panel, a Wristwatch receiver and an LED Receiver and the smaller Cancel Button will cancel the call on multiple receivers.

BENEFITS

SAVE COSTS WITH THIS INTEGRATED BUTTON:

The call and cancel are both treated as one button so if your display has a capacity of 500 then you can have 500 residents with this button by their bedsides. The integrated cancel button will also monitor the response time of the nurse when making use of the reporting software.

PREVENT MULTIPLE CALLS:

The smaller cancel button will light up when the large call button is pressed as confirmation that the call has gone through. This gives the patient or resident peace of mind that they will be attended to and prevents them from pressing the button multiple times in order to 'make sure' that their call went through.

WATERPROOF:

This call Button can be mounted in high-risk areas such as bathrooms, toilets and showers.

Your patient or resident will be able to call out for help immediately with a simple push of a button.

FEATURES

INTEGRATED CANCEL BUTTON:

The integrated cancel button allows the nurse to cancel the call at the call point, whether it is in the patient or resident’s room, toilet or shower.


MOUNTABLE:

This button can be wall mounted with screws and is supplied with a built-in mounting bracket. Alternatively, you can use strong double-sided tape.

BATTERIES:

A23 Battery that can last 1-2 years depending on use.

The LED will flicker at a slower pace and become dull when the battery is low and needs replacement.



You will require a key or flat device to remove the backing.

Insert the key in the smaller of the two indentations to pop off the backing.

Once the backing is removed you can insert the battery, protected by a waterproof rubber flap.

B-SW2 Technical Information	Dimensions: 85x85x15mm	Range: 100 - 150m depending on the surroundings
	Battery: A23	Radio Frequency: 433MHz


CALL AND CANCEL WALL UNIT WITH CORD

A simple push of the button will send out alerts to multiple receivers such as:

- A Display Panel,
- A Wristwatch receiver and
- An LED Receiver.

Can be mounted at the bedside, for frail patients in a Retirement Home or Hospital.

- Comes with an extended Call Button at the end of the cord for easy reach and a more convenient press of the button.
- Can be mounted at the bedside, shower or toilet in a Retirement home or Hospital.
- The large Call Button can send out alerts to multiple receivers such as a Display Panel, a Wristwatch Receiver and an LED Receiver and the smaller Cancel Button will cancel the call on multiple receivers.



BENEFITS

EXTENDED CALL BUTTON:

The call button at the end of the cord can be placed by the bedside for a more convenient press of the button. This unique Call Button has a built-in call and cancel button with an extended Call Button on a cord.

CORD LENGTH:

The cord length is 60 cm (coiled up) and can extend to 1.5 meters when pulled to be used by the bedside.

SAVE COSTS WITH THIS INTEGRATED BUTTON:

The call and cancel are both treated as one button so if your display has a capacity of 500 then you can have 500 residents with this button by their bedsides. The integrated cancel button will also monitor the response time of the nurse when making use of the reporting software.



The B-SW-G33 is mounted flush against the wall or any flat surface with a mounting bracket for the extended button.



The extended button can be placed in the mounting bracket or left on the bed of frail patients.

FEATURES

RESPONSE TIME TRACKING:

This Call button has an integrated cancel button to monitor the Nurse's Response Time when making use of the reporting software.

MOUNTABLE:

This button can be wall mounted with screws and is supplied with a built-in mounting bracket. Alternatively, you can use strong double-sided tape.

BATTERIES:

A23 Battery that can last 1-2 years depending on use.

The LED will flicker at a slower pace and become dull when the battery is low and needs replacement.

B-SW-G33 Technical Information	Dimensions: 85x85x15mm (panel)	Radio Frequency: 433MHz
	Battery: A23	Range: 100 - 150m depending on the surroundings

LED LIGHT RECEIVER

This receiver can be installed above the ward or room entrance situated in the corridors.

It is also the ideal receiver to use as part of a panic button system.



BENEFITS

INSTANT ALERTS:

- Improve your response time to call for help.
- Helps the Nurse to quickly identify the room in need without having to return to the nurse station.

ALERT TONE:

Audio:

The LED has a switch inside the back cover which allows you to change between a continuous beeping noise or silence.

Visual:

This LED light also flashes bright red upon pressing the call button.



LED Light installed above a door frame for easy identification in a corridor.

FEATURES

POWER:

Even though it is part of our wireless range, it does require constant power. The signal transmits wirelessly from the Call Button to the LED light. The LED light is supplied with the required 12V power supply which is 1m long.

B-700

CAPACITY:

This LED Light Receiver can store up to 99 call buttons & requires 1 cancel button. It can be configured by the following two methods:

- 49 Call Buttons and 49 Cancel Buttons
- 98 Call Buttons and 1 Cancel Button



CANCEL WITH ONE OR MULTIPLE:

The LED light cannot work without a Cancel Button.

A cancel button is used to stop the light from flashing and making a noise. You can share one cancel button for multiple call buttons or have a cancel button for each call. The configuration would be dependent on what you like displayed on your display panel, so it is compatible with the display information required.

B-700 Technical Information	Working Voltage: 12V 500mA Working Current: <65mA	Sensitivity: 106dB Volume Output: 85dB	Standby Current: <10mA Working Environment: -40~80°C
	Dimensions: 104x86x38mm	Range: 100 - 150m depending on the surroundings Radio Frequency: 433MHz	

WRISTWATCH RECEIVER

Do you need to receive notifications on the move?

- In a situation where you would like faster response time, not having to go to the Nurse Station to view the details saves you time.
- The wristwatch gives you an immediate notification while attending to another patient or resident.
- The Nurse can cancel the request on the wristwatch, attend to the resident and then cancel at the call point.



FEATURES

WEARABLE ON A LANYARD OR WRIST:

You have two options:

Remove the watch straps and wear the wristwatch on a lanyard or wear it on your wrist.

CAPACITY:

The wristwatch receiver has a capacity of 200 transmitters.

The Buttons can be all Call Buttons and it has a dedicated cancel* button on the wristwatch.

**Note: Pressing the cancel button on the wristwatch will not cancel the call on the display panel. It will only clear the notification on the wristwatch display.*

DISPLAY INFORMATION:

Display either four numbers "0000-9999" or use the first digit as a letter, from A-F. Eg: "A001"- "F999".

RECHARGEABLE BATTERY:

Recharge the lithium-ion battery every 2-5 days in a nursing environment (depending on use). Battery life is indicated on the display.

SIGNAL STRENGTH INDICATION:

You can see the strength of the signal displayed on the wristwatch. If you require it to travel further, you can add signal repeaters for an extended range.

PLEASE NOTE

SINGLE KEY CANCEL BUTTON:

If you use this Cancel Button (B-A1-BC) to attend to the patient or resident you will also have to cancel the call on the wristwatch using its dedicated cancel button.



INTEGRATED CALL & CANCEL BUTTONS:

If you are using the integrated Call and Cancel Buttons (B-SW2 / B-SW-G33), the integration cancel button will also cancel the call on the wristwatch.



Micro USB Charging Port

The wristwatch receiver is charged using the micro-USB charger that is supplied with each watch.

B-650 Technical Information	Working Voltage: 3.7V (Li-ion battery) Working Current: <900mA	Sensitivity: 85dB	Radio Frequency: 433MHz
	Dimensions: 104x86x38mm	Range: 100 - 150m depending on the surroundings	

DISPLAY PANEL

Suitable for Clinics, Retirement homes or specific wards within a Hospital who would like to easily adjust the calling information and additional call buttons as they go.



BENEFITS

CONTROL:

- Control your Display Panel with 4 easy-to-use buttons.
- The B-99P is popular due to the ease of use and adjusting the programming on the front of the Display Panel in 3 easy steps.
- Receive alerts on a large 2-digit display which is easy to read from a distance or with a loud alert when you are not facing the display panel.

FEATURES

DISPLAY OPTIONS:

Simple 2-digit codes to indicate the area the caregiver needs to attend. The B-99P has a 2-digit display. It can display 01-99.

STAND OR MOUNTABLE:

This display panel comes with a tabletop stand.

You can easily place it on your counter, or if you desire it to be mounted it also comes with a mounting bracket.



Stand supplied for table or counter-top display.



Mounting bracket supplied to mount on a wall.

CAPACITY:

Has a capacity of 99 transmitters.

This means that you can store up to 99 Call or Cancel buttons on this display panel.

ALERT TONE:

One distinct ringtone.

The volume level can be adjusted to your preference.

AUTO-CANCEL:

Convenient and easy to cancel.

You can choose between 5-90 seconds, how long you want the calling information to stay on the screen before it cancels automatically.

POWER SUPPLY:

The B-99P is supplied with a DC12V adapter.

LOADSHEDDING:

In the event of loadshedding, a UPS Backup Power Supply can keep the Display Panel powered on.

B-99P Technical Information	Working Voltage: 12V	Working Current: <450mA	Output: 2w Sensitivity: 110dB
	Dimensions: 248x163x23mm	Range: 100 - 150m depending on the surroundings Radio Frequency: 433MHz	

DISPLAY PANEL

Suitable for Clinics, Retirement homes or specific wards within a Hospital who would like to easily adjust the calling information and additional call buttons as they go.



BENEFITS

CONTROL:

- Control your Display Panel with 4 easy-to-use buttons.
- The B-99E is popular due to the ease of use and adjusting the programming on the side of the Display Panel in 3 easy steps.

FEATURES

DISPLAY OPTIONS:

A wide range of codes is available to represent your call areas.

The B-99E has a three-digit display. It can display all numbers, and the first digit can either be the following ten alphabets (A, b, C, d, E, F, H, L, P, U) or numbers only.

EG: A01-U99 and/or 001-999.

ROLLING DISPLAY:

View all calls at a glance.

If multiple calls come through at the same time, it will roll between each call at the desired interval of 5-99 seconds.



ADJUSTABLE VOLUME AND RINGTONES:

Allows you to adjust the volume to suit your environment.

Choose from a range of 7 different ringtone options that best suits your establishment.

CAPACITY:

Each display panel has a capacity of 256 buttons.

The Buttons can be all Call Buttons or a combination of Call and Cancel Buttons depending on how you would like your Nurse Call System to be set up.

QUEUING MODE:

Attend to multiple calls one at a time.

Once you cancel the current call it will show you the next call. You can see how many calls are waiting assistance - the first digit will display the number of outstanding calls.

AUXILIARY OUT:

Need to hear the Display from a distance?

Connect an external HI-FI with an Aux Cable to wake up the neighbours.

CANCEL BUTTON COMPATIBLE:

Choose how your Nurse Call System operates. If you don't want to use the auto-cancel feature the Display Panel supports 2 types of cancels:

- The nurse would have to actively cancel the call button on the screen.
- The nurse will have to go to the call point and cancel the call via a cancel button at the bedside / call point.

AUTO-CANCEL:

Convenient and easy to programme.

You can choose between 5-90 seconds, how long you want the calling information to stay on the screen before it cancels automatically.

POWER SUPPLY:

The B99-E requires a power source and is supplied with a DC12V adapter which is 1m long.



B-99E Technical Information	Working Voltage: 12V	Working Current: <900mA	Output: 2w Sensitivity: 106dB
	Dimensions: 290x158x40mm	Range: 100 - 150m depending on the surroundings Radio Frequency: 433MHz	

DISPLAY PANEL

Suitable for Clinics, Retirement homes or specific wards within a Hospital who would like to easily adjust the calling information and additional call buttons as they go.



BENEFITS

CONTROL:

- Control your Display Panel with 4 easy-to-use buttons.
- The B-99G is popular due to the ease of use and adjusting the programming on the side of the Display Panel in 3 easy steps.

FEATURES

DISPLAY OPTIONS:

A wide range of codes is available to represent your call areas.

The B-99G has a three-digit display. It can display all numbers, and the first digit can either be the following ten alphabets (A, b, C, d, E, F) or numbers only. EG: A01-F99 and/or 001-999.

ROLLING DISPLAY:

View all calls at a glance.

If multiple calls come through at the same time, it will roll between each call at the desired interval of 5-99 seconds.

ADJUSTABLE VOLUME AND RINGTONES:

Allows you to adjust the volume to suit your environment.

Choose from a range of 6 different ringtone options that best suits your establishment.

CAPACITY:

Each display panel has a capacity of 256 buttons.

The Buttons can be all Call Buttons or a combination of Call and Cancel Buttons depending on how you would like your Nurse Call System to be set up.

STANDBY MODE:

When the display panel is on standby mode and no calls are being shown, it will display the date and time.

CANCEL BUTTON COMPATIBLE:

Choose how your Nurse Call System operates. If you don't want to use the auto-cancel feature the Display Panel supports 2 types of cancels:

- The nurse would have to actively cancel the call button on the screen.
- The nurse will have to go to the call point and cancel the call via a cancel button at the bedside / call point.

AUTO-CANCEL:

Convenient and easy to programme.

You can choose between 5-90 seconds, how long you want the calling information to stay on the screen before it cancels automatically.

POWER SUPPLY:

The B99-G requires a power source and is supplied with a DC12V adapter which is 1m long.

B-99G Technical Information	Working Voltage: 12V	Working Current: <900mA	Output: 2w Sensitivity: 106dB
	Dimensions: 232x14x22mm	Range: 100 - 150m depending on the surroundings Radio Frequency: 433MHz	

DISPLAY PANEL

Suitable for Clinics, Retirement homes or Hospitals that have multiple areas or wards to receive alerts from.

- Receive alerts from any of the Call buttons when pressed.
- Record and extract reports of when buttons were pressed and how long it took staff to respond when used with Desktop Software.
- Show detailed information of where a call originates from with 16 characters to display.



Bed 1 CALL
15:05:05

BENEFITS

DISPLAY OPTIONS:

Up to 16 characters that can be numerical or alphabetical.

With this display panel, you can customize the calling information to show the specific ward or bed number or a resident name, etc.

CONTINUOUS DISPLAY:

You have the option to have the calling information displayed continuously until the nurse has pressed the cancel button by the patient or resident bed. This ensures that no calls are missed or forgotten.

FEATURES

ADJUSTABLE VOLUME AND RINGTONES:

Allows you to adjust the volume to suit your environment.

Choose from a range of different ringtone options that best suits your establishment.

CAPACITY:

An unlimited capacity allows you to link as many buttons as required to this display panel via the Desktop Software.

STANDBY MODE:

When the display panel is on standby mode and no calls are being shown, it will display the date and time.

CANCEL BUTTON COMPATIBLE:

Choose how your Nurse Call System operates. If you don't want to use the auto-cancel feature the Display Panel supports 2 types of cancels:

- The nurse would have to actively cancel the call via the Desktop Software.
- The nurse will have to go to the call point and cancel the call via a cancel button at the bedside / call point.

POWER SUPPLY:

The B-901 Display Panel is equipped with an AC100-240V 50/60Hz power cable. This is plugged into a 2-pin adapter and requires constant power.

CONNECTED VIA LAN:

The B-901 Display Panel requires multiple components (of which all need to be connected with LAN Cables):


- The Display
- The Gateway
- A router with an internet connection
- A laptop/PC with the Software

B-901 Technical Information	Working Voltage: 220V	Radio Frequency: 433MHz
	Dimensions: 826x182x44mm	Range: 100 - 150m depending on the surroundings

DESKTOP SOFTWARE

Suitable for Clinics, Retirement homes or Hospitals that have multiple areas or wards to receive alerts from.

- Receive alerts from any of the Call buttons when pressed.
- Record and extract reports of when buttons were pressed and how long it took staff to respond.
- Show detailed information of where a call originates from with 16 characters to display.



BENEFITS

DISPLAY OPTIONS:

Up to 16 characters that can be numerical or alphabetical.

With the software, you can customize the calling information to show the specific ward or bed number or a resident name, etc.

CONTINUOUS DISPLAY:


The call will remain highlighted until it has been cancelled. This ensures that no calls are missed or forgotten.

NO DISPLAY PANEL REQUIRED:


When making use of the W900 Software, you can choose to not have a display panel as the calling information will show on the PC or Laptop when the Software is open.



Software indicating the active calls by highlighted colour.



Show various records such as devices, calls, service types, etc.



View statistics on number of calls, areas, response times, etc.

FEATURES

VOICE ALERT:

When a call button is pressed, you can have the software set to call out the calling information so that if you are not in front of the Laptop or PC, you will still be notified of the alert.

CAPACITY:


An unlimited capacity allows you to link as many buttons as required to this display panel via the Desktop Software.

ON-SCREEN CANCEL

The nurse can cancel a call via the software if you are not making use of cancel buttons at your establishment.

DEVICE MANAGEMENT

Easily discover, authenticate, add, bind, modify or delete Call and/or Cancel buttons via the software.



Easily discover, authenticate, add, bind, modify or delete Call and/or Cancel buttons via the software.



The nurse can cancel a call via the software if you are not making use of cancel buttons at your establishment.

SMART GATEWAY

Suitable for Clinics, Retirement homes or Hospitals that have multiple areas or wards to receive alerts from and used with the W900 Desktop Software.



SIGNAL REPEATER

Your retirement home or hospital may need signal repeaters to carry the signal further.

These signal repeaters are installed where the signal falls short. BellMe has successfully installed up to 11 repeaters on the sites that are 600 x 600m in retirement homes & hospitals with 12 repeaters over three floors spanning 250 meters.



BENEFITS

INDEPENDANCE:

- Give your patient or residents the freedom to move.
- The signal booster extends the range of the Call Button transmission.
- Patients or Residents have complete peace of mind knowing that their Call Button will reach their caregiver whether they are in the Garden, the TV area or the Canteen.

FEATURES

UNLIMITED CAPACITY:

The call buttons do not need to be programmed to the repeater. and thus allows for an unlimited amount of buttons to work with one repeater.

ANTENNAE:

With two Antennae to help receive and transmit the signal, the device ensures that BellMe sufficiently covers your site with ease.

RANGE:

Our systems have a range of about 100 - 150 meters, depending on the surroundings and can be extended with the signal repeater.

POWER SUPPLY:

The Signal repeater is supplied with the required 12V adapter which is 1m long.

POWER:

Even though it is part of our wireless range, all receivers require constant power except the wristwatch receiver.

BATTERY BACKUP:

In case of a power failure, the repeaters have an internal rechargeable battery pack that will power the unit for up to 1.5 - 2 hours.

PLEASE NOTE

HOW DO YOU DETERMINE IF A SIGNAL REPEATER IS NEEDED, WHERE TO PLACE THEM, AND HOW MANY:

Signal repeaters are strategically placed where the signal falls short and each establishment's environment is different depending on various factors. There is no way to estimate exactly how many signal repeaters your establishment will require - this can only be determined when on-site.



Placing a signal repeater where the signal falls short will increase the range further.

EXTENDING THE POWER SUPPLY:

Should your repeater need to be placed in a strategic position which is out of reach of any power points your electrician can extend the Power Supply for you.

B-Q4 Technical Information	Radio Frequency: 433MHz	Range: 100 - 150m depending on the surroundings
----------------------------------	----------------------------	--

USB RECEIVER

The USB Receiver is an optional extra when purchasing any of the smaller display panels. It can also be used in the absence of any display panel. Mainly used to capture the response time of the nurse using the BellMe Reporting Software.



REQUIREMENTS:

- In order to capture the response time with the USB Receiver, you will require any of the BellMe nurse call buttons that have an integrated cancel button.
- When making use of single-key cancel buttons, it will only record when a call button was pressed but not when it was cancelled.
- Windows XP or higher PC/Laptop.


CAPACITY:

This receiver can store up to 110 transmitters (call buttons)


USB INPUT:

This device connects to your PC/Laptop via USB and does not require an external power source.

REPORTING SOFTWARE



Screenshot of the software showing admin login screen which can stay with a Matron or Manager.



Call Date	Call Time	Call Zone	Call Position	Service type	Service personnel	End Time	Average Time	Service state	Incharge Person
06/03/2019	20:04:46	WARD-8-RM2-CALL(0000)	WARD-8-RM2-CALL(0000)	CALL		20:04:08	00:00:17	No answer	
06/03/2019	20:04:51	WARD-8-RM2-CALL(0000)	WARD-8-RM2-CALL(0000)	CALL		20:04:08	00:00:17	1st Warning	
06/03/2019	20:04:19	TV-ROOM-CALL(0001)	TV-ROOM-CALL(0001)	CALL		20:04:08	00:00:44	2nd Warning	
06/03/2019	20:04:24	F-TOILET-1-CALL(0001)	F-TOILET-1-CALL(0001)	CALL		20:04:08	00:00:09	1st Warning	
06/03/2019	20:03:59	F-TOILET-1-CALL(0001)	F-TOILET-1-CALL(0001)	CALL		20:04:08	00:00:09	Attended	
06/03/2019	20:03:36	F-TOILET-1-CALL(0001)	F-TOILET-1-CALL(0001)	CALL		20:04:08	00:00:31	Attended	

Screenshot of the software showing recorded proof of Calling and Cancelling information between patient or resident and caregiver.

FEATURES


ACCOUNTABILITY:

- Depending on your display panel it can be customized so that the nurse on duty cannot cancel a call from the nurse's station or software screen. The nurse will have to cancel from the call point.
- Captures and provides you with recorded proof of the time taken between the patient or resident calling and being attended to.
- Requires administrative logins which will stay with the Matron or Manager so that the staff cannot change the settings on the software.

DETAILED REPORTS:

When viewing reports you will have access to:

- When a call button was pressed
- When the cancel button was pressed
- Floor plan of your layout (if set up)
- Detailed calling information
- Statistics and graphs on staff response time
- Export reports to Excel for record-keeping



A	B	C	D	E	F	G	H
Date	Call Time	Call Position	Service type	Average Time	Service state	End Time	
9/08/03	20:03:36	F-TOILET-1-CALL(0001)	CALL	0:00:31	Attended	20:04:08	
9/08/03	20:03:59	F-TOILET-1-CALL(0001)	CALL	0:00:09	Attended	20:04:08	
9/08/03	20:04:09	TV-ROOM-CALL(0003)	CALL	0:00:50	Final Warning	0:00:00	
9/08/03	20:04:19	TV-ROOM-CALL(0003)	CALL	0:00:44	2nd Warning	0:00:00	
9/08/03	20:04:31	WRD-C RM10-CALL(0008)	CALL	0:00:32	1st Warning	0:00:00	
9/08/03	20:04:46	WARD-B RM2-CALL(0006)	CALL	0:00:17	No Answer	0:00:00	

Export your reports to excel for record-keeping.

CLIENT

TESTIMONIALS



Christo J Hattingh (Manager)
KarooSIG Aftreeoord - Orania

We purchased this product because of its good price and simple operation. We specifically like that this product has good range, is expandable and upgradable.

The product was affordable and we had quick service.



Bongi Zuma (Acting Study Coordinator: HVTN)
CAPRISA eThekweni Research Centre - Durban


The Nurse Call System assisted us in speeding up the processes in our Clinic setting. If people have buzzed for assistance it displays in the common room. It is good value for money especially if you want to streamline processes and procedures in the clinic.

We were very happy with the excellent service and patience of the staff.



Vicki Bannister (Manager)
Clarendon House - Pietermaritzburg

The product we liked the most was the wristwatch. The system is very user and staff-user friendly and the BellMe staff are very helpful.



Wilma Smit
Lalamanzi - Middelburg

What we liked the most was the fact that we could even pick up bells in the 300m distance. The best product ever. The after sales is excellent and delivery on time.

CLIENT

TESTIMONIALS ”



Justus Vermaak (Director)
Henley River Care Centre - Alberton

Our primary reason for buying this system was for use by frail-care residents to alert nursing staff. We are very happy with the reliability, range of signal and ease of use. It's definitely a good product and an excellent way to alert staff to any arising situation. We were very happy with BellMe's excellent sales and after sales service received.



Sr. Olive Kumalo (Registered Nurse)
Oasis Recovery Centre - Plettenberg Bay

The Nurse Call System provides communication between our patients (from their rooms and bedrooms) and our nurses to ensure our client's well-being at all times. We have peace of mind knowing that we are always available to our clients and knowing exactly who needs to help where. You can't go wrong with this system. It works really well and we have had no problems with it.



Petro Fourie (Manager)
Kenwyn Retirement Home - Pietermaritzburg

Ushir was great in assisting us with all the programming of the buttons. All my residents are now wearing their buttons around their necks...Before handing the buttons out we had a few tested from different rooms and they all came through without any problems. We are very happy with the system so far and would like to thank all involved for the speedy and efficient way things were dealt with.

CLIENT

TESTIMONIALS ”



Sarah Clark (Facilities Manager)
Queenshaven Village - Johannesburg

It has proven to be very useful and as it is extremely simple to use our residents do not have to learn complicated systems which they might struggle to remember in an emergency. It has been put to very effective use on a great number of occasions providing quick responses to calls for help.



Matron
Serene Park Retirement Centre - Pretoria




















The system is very very very helpful because when it rings you can see the letters very well, they are big and clear and the sound is very very clear. Even if you are far, you can hear it and it rings until you come and stop it and then you go to the resident and help them. There are no problems at all because there is no way that you can miss it.



Sr Leslie Bossert (Nurse Practitioner)
Villa Sunfield - Durban

BellMe never failed to deliver. The service I received was excellent and the product is of very high quality. Our residents now are in safer hands knowing that emergency assistance is available at the touch of a button.

SOME OF OUR CLIENTS ENJOYING THE BELLME NURSE CALL SYSTEM

		
		
		
		
		
		
	<p>AKESO BRANCHES SUPPLIED:</p> <ul style="list-style-type: none"> • Randburg – Crescent Clinic • Pietermaritzburg • Kenilworth Adolescents • Umhlanga • Kenilworth Clinic • Stepping Stones • Milnerton • Nelspruit 	

THANK YOU



Please feel free to get in touch with us for any and all enquiries. We would love to hear from you.

AGENTS CURRENTLY AVAILABLE IN THESE AREAS:

- Kwa-Zulu Natal
- Western Cape
- Free State
- North West
- Northern Cape
- Eastern Cape
- Gauteng
- Mpumalanga
- Zimbabwe

DIRECTOR
Ushir Shah



MANAGER
Jolandie van Zyl



CONTACT US



PHONE

+27 64 532 0828



EMAIL

info@bellme.co.za



WEBSITE

www.bellme.co.za